

# Analysis of urine samples using the MxP<sup>®</sup> Quant 500 kit

## Content

1	Introduction and additional items.....	1
2	MetIDQ™ software – Generate worklist in MetLIMS .....	3
3	Instrumental setup.....	4
4	Preparing the kit in the lab.....	4
5	MetIDQ™ software – Normalization and data export .....	6
6	Analytical specifications for creatinine .....	6

## 1 Introduction and additional items

A protocol for the analysis of urine samples using the MxP<sup>®</sup> Quant 500 kit is described below. The workflow is slightly different from the one described in the user manual and additional items are required, as listed in the table below.

In general, we recommend performing pilot tests with a representative sample set before starting a larger study. The results may vary depending on the nature and the quality of the samples. Due to the high salt content of urine samples, it is recommended to clean the instrument after every kit. biocrates does not assume responsibility for the results or possible system contamination.

The following items are required in addition to the standard MxP® Quant 500 kit contents

Additional items	Description	Details
<b>Triethylamine</b> <b>Purchased by customer</b>	Purity ≥99%	Used for derivatization solution <b>instead</b> of pyridine
<b>MetIDQ™ patch:</b> MStype_Q500_urine_DB110_YYMMDD.jar (provided by support)	Patch containing urine analysis protocols for MetIDQ™	To be imported into MetIDQ™ <b>in addition</b> to the regular Quant 500 patch, required to generate the worklist
<b>AbsoluteIDQ® p180 Cal1 – Cal7</b> , 7 vials	Calibration standards, lyophilized	Used <b>instead</b> of the regular Quant 500 calibration standards
<b>AbsoluteIDQ® p150/p180 ISTD Urine Crea</b> , 1 vial	Internal standard creatinine, lyophilized	<b>Additional</b> creatinine internal standard to be added to plate
<b>AbsoluteIDQ® p150/p180 Urine Zero Sample</b> , 1 vial	Solution of 350 mM urea, 15 mM NaHPO <sub>4</sub> , pH 6.0 (urine imitation)	Used as zero sample <b>instead</b> of PBS
<b>Recipe urine control</b> , 1 vial	Lyophilized urine QC sample	Used as <b>additional</b> quality control
<b>Methods for Agilent MassHunter software</b> (provided by support)	<b>LC1:</b> MxP500L-LC1_7x11_urine.m <b>FIA1:</b> MxP500F-FIA1_7x11_urine.m <b>LC2, FIA2 and FIA3:</b> regular according to manual	Used for data acquisition. <b>The specific numbers (x) depend on MS type, see Quant 500 and MetIDQ™ user manuals.</b>
<b>Methods for Waters MassLynx® software</b> (provided by support)	<b>LC1:</b> MxP500L-LC1_8xx2_urine.exp <b>FIA1:</b> MxP500F-FIA1_8xx2_urine.exp <b>LC2 and FIA2 (and FIA3 for TQ-S):</b> regular according to manual	Used for data acquisition. <b>The specific numbers (xx) depend on MS type, see Quant 500 and MetIDQ™ user manuals.</b>
<b>Methods for SCIEX Analyst® software</b> (provided by support)	<b>LC1:</b> MxP500L-LC1_5xx2_urine.dam <b>FIA1:</b> MxP500F-FIA1_5xx2_urine.dam <b>LC2 and FIA2:</b> regular according to manual	Used for data acquisition. <b>The specific numbers (xx) depend on MS type, see Quant 500 and MetIDQ™ user manuals.</b>

## 2 MetIDQ™ software – Generate worklist in MetLIMS

Step	Instruction																																	
1	Install the patch <b>MStype_Q500_urine_DB110_YYMMDD.jar</b> (provided by support) according to MetIDQ™ user manual section “Installing database patches”.																																	
2	Select the material class “(40) urine”, when you register urine samples.																																	
3	Select the correct OP for your instrument when you generate the LC worklist for urine: <b>MXP500L-40-xxxx</b>																																	
4	In the <b>Zero Samples</b> tab of the <b>Worklists Settings For Submission</b> window, unlink PBS and link <b>Urine imitation</b> (barcode 11000007).																																	
5	In the <b>Standards</b> tab of the <b>Worklists Settings For Submission</b> window, link the following calibration standards: <table border="1" data-bbox="335 884 1273 1444"> <thead> <tr> <th>Sample Bar Code</th> <th>Sample Type</th> <th>Sample Identification</th> </tr> </thead> <tbody> <tr> <td>1020509951</td> <td>Standard L0.25</td> <td>p180 Cal0.25</td> </tr> <tr> <td>1020509966</td> <td>Standard L0.5</td> <td>p180 Cal0.5</td> </tr> <tr> <td>20000611</td> <td>Standard L1</td> <td>p180 Cal1</td> </tr> <tr> <td>20000612</td> <td>Standard L2</td> <td>p180 Cal2</td> </tr> <tr> <td>20000613</td> <td>Standard L3</td> <td>p180 Cal3</td> </tr> <tr> <td>20000614</td> <td>Standard L4</td> <td>p180 Cal4</td> </tr> <tr> <td>20000615</td> <td>Standard L5</td> <td>p180 Cal5</td> </tr> <tr> <td>20000616</td> <td>Standard L6</td> <td>p180 Cal6</td> </tr> <tr> <td>20000617</td> <td>Standard L7</td> <td>p180 Cal7</td> </tr> <tr> <td>1020509971</td> <td>Standard L8</td> <td>p180 Cal8</td> </tr> </tbody> </table>	Sample Bar Code	Sample Type	Sample Identification	1020509951	Standard L0.25	p180 Cal0.25	1020509966	Standard L0.5	p180 Cal0.5	20000611	Standard L1	p180 Cal1	20000612	Standard L2	p180 Cal2	20000613	Standard L3	p180 Cal3	20000614	Standard L4	p180 Cal4	20000615	Standard L5	p180 Cal5	20000616	Standard L6	p180 Cal6	20000617	Standard L7	p180 Cal7	1020509971	Standard L8	p180 Cal8
Sample Bar Code	Sample Type	Sample Identification																																
1020509951	Standard L0.25	p180 Cal0.25																																
1020509966	Standard L0.5	p180 Cal0.5																																
20000611	Standard L1	p180 Cal1																																
20000612	Standard L2	p180 Cal2																																
20000613	Standard L3	p180 Cal3																																
20000614	Standard L4	p180 Cal4																																
20000615	Standard L5	p180 Cal5																																
20000616	Standard L6	p180 Cal6																																
20000617	Standard L7	p180 Cal7																																
1020509971	Standard L8	p180 Cal8																																
6	In the <b>QC-Lots</b> tab of the <b>Worklists Settings For Submission</b> window, link the regular <b>Quant 500 QC levels 1-3</b> , as well as the <b>Recipe Urine QC</b> (QC Level 0). Adjust at least <b>3 replicates</b> of both, Quant 500 <b>QC2</b> and Recipe Urine QC, and distribute over the plate.																																	
7	In the <b>Samples</b> tab of the <b>Worklists Settings For Submission</b> window, link the experimental samples as usual.																																	
8	Duplicate the LC plate (copy and paste derived plate) as usual and change the OP of the copied plate to the corresponding one for FIA: <b>MXP500F-40-xxxx</b>																																	
9	According to the regular Quant 500 workflow (see user manual), delete the calibration standards in the <b>Plate View</b> of the FIA plate, set the <b>Condition</b> to <b>Approved</b> and perform <b>Export Worklist for MS</b> .																																	

### 3 Instrumental setup

Step	Instruction
1	On the MS operating computer, copy the instrument-specific acquisition methods (see page 1 of this document) into the corresponding folder of your biocrates project.
2	Install the acquisition methods according to section 4 of the MxP® Quant 500 kit user manual.
3	Continue with the MxP® Quant 500 kit user manual and run the SST and the kit.



Do **not** thaw samples or start with the kit before successful completion of the SST.

### 4 Preparing the kit in the lab


Step	Instruction												
1	Resuspend all vials using LC-MS grade water:												
	<table border="1"> <thead> <tr> <th>Item</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>ISTD Urine Crea</td> <td>add 1200 µL</td> </tr> <tr> <td>AbsolutelDQ® p180 Cal1 - Cal6</td> <td>add 100 µL</td> </tr> <tr> <td>AbsolutelDQ® p180 Cal7</td> <td>add 50 µL</td> </tr> <tr> <td>MxP® Quant 500 QC1 - QC3</td> <td>add 100 µL</td> </tr> <tr> <td>Recipe Urine QC</td> <td>add 1000 µL</td> </tr> </tbody> </table>	Item	Volume	ISTD Urine Crea	add 1200 µL	AbsolutelDQ® p180 Cal1 - Cal6	add 100 µL	AbsolutelDQ® p180 Cal7	add 50 µL	MxP® Quant 500 QC1 - QC3	add 100 µL	Recipe Urine QC	add 1000 µL
	Item	Volume											
	ISTD Urine Crea	add 1200 µL											
	AbsolutelDQ® p180 Cal1 - Cal6	add 100 µL											
	AbsolutelDQ® p180 Cal7	add 50 µL											
MxP® Quant 500 QC1 - QC3	add 100 µL												
Recipe Urine QC	add 1000 µL												
2	Shake all vials for 15 min at 1200 rpm and vortex several times before use.												
3	Add 10 µL of the <b>ISTD Urine Crea</b> to each well of the kit plate <b>except the blank well A1</b> . Pipette directly onto the filters of the kit plate. Do not pipette on the inner wall of the wells or on the plastic holder. Use a repeater, e.g. Eppendorf Multipipette®, adjusted to maximum dispensing speed.												

Step	Instruction																						
4	<p>Load the resuspended calibration standards from step 1 directly on the kit plate as follows (according to your plate map):</p> <table border="1" data-bbox="336 504 1023 1059"> <thead> <tr> <th>Calibration level</th> <th>Volume from vial "p180 Cal"</th> </tr> </thead> <tbody> <tr> <td>Cal 0.25</td> <td>2.5 <math>\mu</math>L of p180 Cal 1</td> </tr> <tr> <td>Cal 0.5</td> <td>5 <math>\mu</math>L of p180 Cal 1</td> </tr> <tr> <td>Cal 1</td> <td>10 <math>\mu</math>L of p180 Cal 1</td> </tr> <tr> <td>Cal 2</td> <td>10 <math>\mu</math>L of p180 Cal 2</td> </tr> <tr> <td>Cal 3</td> <td>10 <math>\mu</math>L of p180 Cal 3</td> </tr> <tr> <td>Cal 4</td> <td>10 <math>\mu</math>L of p180 Cal 4</td> </tr> <tr> <td>Cal 5</td> <td>10 <math>\mu</math>L of p180 Cal 5</td> </tr> <tr> <td>Cal 6</td> <td>10 <math>\mu</math>L of p180 Cal 6</td> </tr> <tr> <td>Cal 7</td> <td>5 <math>\mu</math>L of p180 Cal 7</td> </tr> <tr> <td>Cal 8</td> <td>10 <math>\mu</math>L of p180 Cal 7</td> </tr> </tbody> </table> <p>Use a single-channel pipette to pipette the volumes according to the table directly onto the center of each filter. Gently touch the filter inserts with the pipette tip while pipetting the samples. Do not pipette on the inner wall of the wells or on the plastic holder and avoid cross-contamination. Use a fresh tip for each sample.</p>	Calibration level	Volume from vial "p180 Cal"	Cal 0.25	2.5 $\mu$ L of p180 Cal 1	Cal 0.5	5 $\mu$ L of p180 Cal 1	Cal 1	10 $\mu$ L of p180 Cal 1	Cal 2	10 $\mu$ L of p180 Cal 2	Cal 3	10 $\mu$ L of p180 Cal 3	Cal 4	10 $\mu$ L of p180 Cal 4	Cal 5	10 $\mu$ L of p180 Cal 5	Cal 6	10 $\mu$ L of p180 Cal 6	Cal 7	5 $\mu$ L of p180 Cal 7	Cal 8	10 $\mu$ L of p180 Cal 7
Calibration level	Volume from vial "p180 Cal"																						
Cal 0.25	2.5 $\mu$ L of p180 Cal 1																						
Cal 0.5	5 $\mu$ L of p180 Cal 1																						
Cal 1	10 $\mu$ L of p180 Cal 1																						
Cal 2	10 $\mu$ L of p180 Cal 2																						
Cal 3	10 $\mu$ L of p180 Cal 3																						
Cal 4	10 $\mu$ L of p180 Cal 4																						
Cal 5	10 $\mu$ L of p180 Cal 5																						
Cal 6	10 $\mu$ L of p180 Cal 6																						
Cal 7	5 $\mu$ L of p180 Cal 7																						
Cal 8	10 $\mu$ L of p180 Cal 7																						
5	<p>Load <b>10 <math>\mu</math>L</b> of all other samples (zero sample, QCs and study samples) as usual according to the user manual and your plate map. Use a single-channel pipette to pipette 10 <math>\mu</math>L onto the center of each filter. Gently touch the filter inserts with the pipette tip while pipetting the samples. Do not pipette on the inner wall of the wells or on the plastic holder and avoid cross-contamination. Use a fresh tip for each sample.</p>																						
6	<p>Dry all wells for <b>30 min</b> under nitrogen according to the user manual.</p>																						
7	<p><b>Derivatization (<u>different</u> from the regular manual):</b>          Prepare the derivatization pre-mix by adding the following chemicals to the plastic tube included in the kit box:</p> <ul style="list-style-type: none"> <li>– 4.2 mL of methanol (LC-MS grade)</li> <li>– 0.6 mL of water (LC-MS grade)</li> <li>– 0.6 mL of triethylamine (<math>\geq</math>99% purity)</li> </ul> <p>Vortex for 10 sec.</p>																						
8	<p>Remove the phenyl isothiocyanate (PITC) from the freezer and allow to equilibrate to room temperature. Prepare the derivatization solution by adding <b>0.6 mL</b> of PITC to the derivatization pre-mix. Vortex for 10 sec.</p>																						
9	<p>Add <b>50 <math>\mu</math>L</b> of the derivatization solution to each well using a repeater, e.g. Eppendorf Multipette®. The derivatization time at room temperature is <b>20 min</b>.</p>																						
10	<p>Continue with the standard manual protocol (Step 7 – Quant 500 manual) and dry all wells for <b>60 min</b> under nitrogen.</p>																						

## 5 MetIDQ™ software – Normalization and data export

Data normalization can be performed using the **Recipe Urine QC**. Furthermore, creatinine is used specifically for normalization of urine concentrations (please refer to Waikar et al., *Kidney Int* 2010; 78(5):486-94). Both normalizations can be automatically performed in MetIDQ™.

Step	Instruction
1	Go to <b>MetSTAT &gt; Select Samples</b> and link the samples of all plate runs of your urine project. Make sure there are two lines for each linked sample with plate barcode and OP belonging to the LC and FIA injection, respectively (except calibration standards which are only measured in the LC part). This will merge the LC and FIA data what is <b>required</b> for the creatinine normalization (step 3).
2	Go to the <b>Display Data</b> tab and use the normalization settings in the tool <b>Data Normalization</b> on the right sidebar as shown in the screenshot below. In this case, QC Level 0 is the <b>Recipe Urine QC</b> . <b>Important:</b> Do not normalize using QC Level 2, which is the regular Quant 500 QC2 (plasma-based), as this will significantly diminish the urine results.
3	Activate the checkbox <b>Creatinine Normalization</b> and each metabolite concentration will be automatically divided by the creatinine concentration.



**Data Normalization**

Normalize Sample Concentration Data Show Report

Plate Source: All

Sample Source: QC Level 0

Normalization vs. Target Values:

Method: Mean

Log Transform Data log<sub>2</sub>

**Creatinine Normalization**

## 6 Analytical specifications for creatinine

Analyte	LOD (µM)	LLOQ (µM)	ULOQ (µM)
Creatinine	100	500	30,000



Please check our [support FAQ](#) for solutions to common performance issues and technical problems or contact us: [support@biocrates.com](mailto:support@biocrates.com).