


Gel Band Submission Form



Name: _____ UM X.500 _____

Date: _____

Species: _____

Gel type: _____

Protein-staining reagent: _____

Amount of protein loaded into each lane (user additional pages if necessary):

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

Lane _____ ug protein loaded _____

MW markers/protein standard (label on gel and supply name/type of marker here):

REQUIRED: Submit gel image with labels as a file or printed copy

see Example "how to label your gel image" here: <http://cbs.umn.edu/cmsp/sample-submission>

Samples: On gel image, gel regions marked in boxes and labeled with sample names were excised and submitted for proteolytic digestion and LC-MS analyses

Enzyme: trypsin _____ or other _____